

## **Genetically modified DBN-09004-6 Soya bean seed powder (blank)**

Art. ID                      ERM-BF441a  
Unit                         vial  
Deliverydetails            No Dangerous Good /not restricted

### Description

ANALYTICAL METHODS USED FOR CHARACTERISATION: Gravimetry. Event-specific quantitative polymerase chain reaction (qPCR).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Genetically modified soya bean with the unique identifier DBN-09004-6	DBN-09004-6		0,00 + 0,08 - 0,00	g/kg		