

Genetically modified DBN-09004-6 Soya bean seed powder (nominal 0.1 %)

Art. ID ERM-BF441c
Unit vial
Deliverydetails No Dangerous Good /not restricted

Description

ANALYTICAL METHODS USED FOR CHARACTERISATION: Gravimetry. Event-specific quantitative polymerase chain reaction (qPCR).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Genetically modified soya bean with the unique identifier DBN-09004-6	DBN-09004-6		1,00 ± 0,11	g/kg		