

### **Custom Solution - 7 components 500 mg/L each of in Methanol**

Art. ID LM24-CUS-72113-3X1ML  
Unit 3 x 1 mL ampoule  
Deliverydetails De Minimis UN1230 (6.1) II E2

#### Description

Conservation 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Benzene	[71-43-2]	500	mg/L		
	Toluene	[108-88-3]	500	mg/L		
	Ethylbenzene	[100-41-4]	500	mg/L		
	o-Xylene	[95-47-6]	500	mg/L		
	m-Xylene	[108-38-3]	500	mg/L		
	p-Xylene	[106-42-3]	500	mg/L		
	Styrene	[100-42-5]	500	mg/L		