

### **Custom Solution - 2 components in HNO3 2%**

Art. ID LM24-CUS-72115  
Unit 5 x 100 mL bottle  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

#### Description

Conservation 24 months. (12 months after opening)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	P from Ammonium Dihydro gen Phosphate		100	mg/L		
	Si from Silicon Dioxide		1000	mg/L		