

Custom Solution - 6 components in HNO3 5%

Art. ID LM24-CUS-72116
Unit 5 x 100 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation 24 months. (12 months after opening)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Ca from Calciumnitrate		5000	mg/L		
	Na from Sodiumnitrate		1000	mg/L		
	Mg from Magnesiumnitrat e		1000	mg/L		
	K from Potassiumnitrate		100	mg/L		
	B from Boricacid		100	mg/L		
	Ba from Bariumnitrate		100	mg/L		