

## **Sulfur in Gasoline High Level Calibration Set**

Art. ID	ACCU-D-7039-HL-SET
Unit	6 x 1 mL
Deliverydetails	De Minimis UN1993 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
ACCU-D-7039-HL-01	Sulfur (S)	[7704-34-9]	100	µg/g		
ACCU-D-7039-HL-02	Sulfur (S)	[7704-34-9]	200	µg/g		
ACCU-D-7039-HL-03	Sulfur (S)	[7704-34-9]	300	µg/g		
ACCU-D-7039-HL-04	Sulfur (S)	[7704-34-9]	400	µg/g		
ACCU-D-7039-HL-05	Sulfur (S)	[7704-34-9]	500	µg/g		
ACCU-D-7039-HL-BL, blank	Sulfur (S)	[7704-34-9]	0	µg/g		

k