

Custom Solution - 5 components in HNO3 5%

Art. ID LM24-CUS-72127
Unit 1 x 100 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation 12 months. (6 months after opening)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Zinc (Zn)	[7440-66-6]	1000	mg/L		
	Iron (Fe)	[7439-89-6]	500	mg/L		
	Lead (Pb)	[7439-92-1]	500	mg/L		
	Vanadium (V)	[7440-62-2]	500	mg/L		
	Silver (Ag)	[7440-22-4]	50	mg/L		