

### **Custom Solution - Phosphates as P (P (PO43-)) 3 mg/L in Water**

Art. ID                      LM24-CUS-72138  
Unit                         5 x 100 mL bottle  
Deliverydetails            No Dangerous Good /not restricted

#### Description

Conservation 12 months. (6 months after opening)

| Text/Information             | Analyte/Parameter           | CAS number | Concentration/Value | Unit | Method | Source |
|------------------------------|-----------------------------|------------|---------------------|------|--------|--------|
| Phosphates as P (P (PO4 3-)) | Phosphates as P (P(PO43 -)) |            | 3                   | mg/L |        |        |