

2',3'-cGAMP, Ammonium Salt (in Solution) (Adenosine-1',2',3',4',5',5"-D6, 98%) (chemical purity 90%)

Art. ID CIL-DLM-11273-CA-25
Unit 25 µmol
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2',3'-cGAMP, Ammonium S alt (Adenosine-1',2',3' ,4',5',5"-D6, 98%)					