

Dabigatran Heptane-2-yl Analog (ethyl 3-(2-(((4-(N-((heptan-2-yloxy)carbonyl)carbamimidoyl)phenyl)amino)methyl)-1-methyl-N-(pyridin-2-yl)-1H-benzo[d]imidazole-5-carboxamido)propanoate)

Art. ID USP-1A03750

Unit 25 mg

Description

Note: PAI products are different from official USP reference standards. These products are released using a separate process developed by USP's subject matter experts. The release process follows rigorous internal policies, standard operating procedures, and USP Quality Management System requirements. USP is an ISO9001:2015 certified facility.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Dabigatran Heptane-2-yl Analog	[1610758-21-8]				