

4-Nitrosopiperazine-1-ethanol Solution (1 mg/mL in Methanol)
(2-(4-Nitrosopiperazin-1-yl)ethan-1-ol)

Art. ID USP-1A04100
Unit 1 mL
Deliverydetails UN1230 (6.1) II E2

Description

Note: PAI products are different from official USP reference standards. These products are released using a separate process developed by USP's subject matter experts. The release process follows rigorous internal policies, standard operating procedures, and USP Quality Management System requirements. USP is an ISO9001:2015 certified facility.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	4-Nitrosopiperazine-1-e thanol	[48121-20-6]	1	mg/mL		