

Phosphatidylethanol 16:0/18:1-D5 (PEth 16:0/18:1-D5) solution 100 µg/mL in methanol ((as free phosphate)), certified reference material, Cerilliant

Art. ID SAF-P-164-1ML
Unit 1 mL
Deliverydetails Dry ice shipment UN1230 (6.1) II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Phosphatidylethanol 16:0/18:1-D5 (PEth 16:0/18:1-D5)	[1612140-92-7]	100	µg/mL		