

Custom Solution - 11 components in HNO3 5%

Art. ID LM24-CUS-72158
Unit 1 x 100 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation 12 months. (6 months after opening)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Bismuth (Bi)	[7440-69-9]	2,5	mg/L		
	Cadmium (Cd)	[7440-43-9]	2,5	mg/L		
	Chromium (Cr)	[7440-47-3]	2,5	mg/L		
	Copper (Cu)	[7440-50-8]	2,5	mg/L		
	Iron (Fe)	[7439-89-6]	2,5	mg/L		
	Magnesium (Mg)	[7439-95-4]	2,5	mg/L		
	Manganese (Mn)	[7439-96-5]	2,5	mg/L		
	Lead (Pb)	[7439-92-1]	2,5	mg/L		
	Barium (Ba)	[7440-39-3]	100	mg/L		
	Calcium (Ca)	[7440-70-2]	100	mg/L		
	Strontium (Sr)	[7440-24-6]	100	mg/L		