

Custom Solution - 9 components 100 mg/L each of in HNO3 10%

Art. ID LM24-CUS-72164
Unit 100 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN2031 II E2

Description

Conservation 24 months. (12 months after opening)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Arsenic (As)	[7440-38-2]	100	mg/L		
	Cadmium (Cd)	[7440-43-9]	100	mg/L		
	Copper (Cu)	[7440-50-8]	100	mg/L		
	Iron (Fe)	[7439-89-6]	100	mg/L		
	Manganese (Mn)	[7439-96-5]	100	mg/L		
	Mercury (Hg)	[7439-97-6]	100	mg/L		
	Lead (Pb)	[7439-92-1]	100	mg/L		
	Tin (Sn)	[7440-31-5]	100	mg/L		
	Zinc (Zn)	[7440-66-6]	100	mg/L		