

### **Custom Solution - 7 components in HNO3 2%**

Art. ID LM24-CUS-72174  
Unit 2 x 250 mL bottle  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

#### **Description**

Conservation 8 months. (4 months after opening)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	50	µg/L		
	Arsenic (As)	[7440-38-2]	10	µg/L		
	Chromium (Cr)	[7440-47-3]	25	µg/L		
	Iron (Fe)	[7439-89-6]	50	µg/L		
	Manganese (Mn)	[7439-96-5]	20	µg/L		
	Nickel (Ni)	[7440-02-0]	20	µg/L		
	Phosphorus (P)	[7723-14-0]	10	µg/L		