

**Liquid Analytical Standard - 100 µg/mL Deoxynivalenol in Methanol**

Art. ID TRIL-TAS-M21LM2-10

Unit 10 mL

Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

| Text/Information | Analyte/Parameter | CAS number   | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|--------------|---------------------|------|--------|--------|
|                  | Deoxynivalenol    | [51481-10-8] | 100                 | µ/mL |        |        |