

Custom Solution - 8 Components 100 mg/L each of in Dichloromethane

Art. ID LM24-CUS-72184
Unit 3 x 1 mL ampoule
Deliverydetails De Minimis UN1593 III E1

Description

Conservation 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	alpha-Chlorotoluene	[100-44-7]	100	mg/L		
	alpha-alpha-Dichlorotoluene	[98-87-3]	100	mg/L		
	alpha,alpha,alpha-Trichlorotoluene	[98-07-7]	100	mg/L		
	alpha,alpha,alpha-4-Tetrachlorotoluene	[5216-25-1]	100	mg/L		
	alpha,alpha,alpha-2-Tetrachlorotoluene	[2136-89-2]	100	mg/L		
	alpha,alpha-2,6-Tetrachlorotoluene	[81-19-6]	100	mg/L		
	3,4-Dichlorobenzotrichloride	[13014-24-9]	100	mg/L		
	2,4-Dichlorobenzyl dichloride	[134-25-8]	100	mg/L		