

EPA Method 533 Native Analyte PDS in methanol (with 4 molar equivalent NaOH)

Art. ID CIL-ES-5654
Unit 1,2 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Perfluorobutyric acid (PFBA)	[375-22-4]	1000	ng/mL		
	Perfluoropentanoic acid (PFPeA), sodium salt	[2706-89-0]	1000	ng/mL		
	Perfluorohexanoic acid (PFHxA), sodium salt	[2923-26-4]	1000	ng/mL		
	Perfluoroheptanoic acid (PFHpA)	[375-85-9]	1000	ng/mL		
	Perfluorooctanoic acid (PFOA)	[335-67-1]	1000	ng/mL		
	Perfluorononanoic acid (PFNA)	[375-95-1]	1000	ng/mL		
	Perfluorodecanoic acid (PFDA), sodium salt	[3830-45-3]	1000	ng/mL		
	Perfluoroundecanoic acid (PFUnA), sodium salt	[60871-96-7]	1000	ng/mL		
	Perfluorodecanoic acid (PFDoA), sodium salt	[3830-45-3]	1000	ng/mL		
	Perfluorobutanesulfonate (PFBS), potassium salt	[29420-49-3]	1000	ng/mL		
	Perfluoropentanesulfonate (PFPS), sodium salt	[630402-22-1]	1000	ng/mL		
	Perfluorohexanesulfonate (PFHxS), potassium salt (mixed isomers)	[3871-99-6]	1000	ng/mL		
	Perfluoroheptanesulfonate (PFHpS), sodium salt	[21934-50-9]	1000	ng/mL		
	Perfluorooctanesulfonate (PFOS) (mixed isomers)	[1763-23-1]	1000	ng/mL		

1H,1H,2H,2H-Perfluorohe xanesulfonate (4:2 FTS) , sodium salt	[27619-93-8]	1000	ng/mL
1H,1H,2H,2H-Perfluoroc tanesulfonate (6:2 FTS) , sodium salt	[27619-94-9]	1000	ng/mL
1H,1H,2H,2H-Perfluorode canesulfonate (8:2 FTS) , sodium salt	[27619-96-1]	1000	ng/mL
Tetrafluoro-2-(heptaflu oropropoxy)propanoic ac id (HFPO-DA) ("GenX")	[13252-13-6]	1000	ng/mL
Dodecafluoro-3H-4,8-dio xanonanoic acid (DONA)	[919005-14-4]	1000	ng/mL
9-Chlorohexadecafluoro- 3-oxanone-1-sulfonic ac id (9Cl-PF3ONS), potass ium salt	[73606-19-6]	1000	ng/mL
11-Chloroeicosafuoro-3 -oxaundecane-1-sulfonic acid (11Cl-PF3OUdS), p otassium salt	[83329-89-9]	1000	ng/mL
Perfluoro-3-methoxyprop anoic acid (PFMPA)	[377-73-1]	1000	ng/mL
Perfluoro-4-methoxybuta noic acid (PFMBA)	[863090-89-5]	1000	ng/mL
Perfluoro-3,6-dioxahept anoic acid (PFDHA)	[151772-58-6]	1000	ng/mL
Perfluoro(2-ethoxyethan e)sulfonic acid (PFEESA)	[113507-82-7]	1000	ng/mL