

Plastic Polyethylene Check Standards for WEEE RoHS Requirements, Disc Ø 31 mm x thickness on request

Art. ID ASI-ROHS-PE9-5ED-QC-31
Unit disc
Deliverydetails No Dangerous Good /not restricted

Description

QC Standard for ASI-ROHS-PE9-5ED-31. Plastic Polyethylene Calibration Standards for WEEE RoHS Requirements, Disc, 5 elements per standard, 31 mm x thickness on request, Br, Cd, Cr, Hg, Pb.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
QC Sample 1	Bromine (Br)	[7726-95-6]	0,025	%		
QC Sample 1	Cadmium (Cd)	[7440-43-9]	0,005	%		
QC Sample 1	Chromium (Cr)	[7440-47-3]	0,05	%		
QC Sample 1	Mercury (Hg)	[7439-97-6]	0,05	%		
QC Sample 1	Lead (Pb)	[7439-92-1]	0,05	%		