

Custom Solution - Si 2000 mg/L in HNO3 5%

Art. ID LM24-CUS-72197
Unit 1 x 100 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation 24 months. (12 months after opening)

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Silicon (Si) | [7440-21-3] | 2000 | mg/L | | |