

Custom Solution - 5 components in 0.5 % Nitric Acid

Art. ID LM24-CUS-72201
Unit 1 x 100 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG)

Description

Conservation 24 months. (12 months after opening)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Sodium (Na ⁺)		2500	mg/L		
	Ammonia (NH ₃)		2500	mg/L		
	Potassium (K ⁺)		1250	mg/L		
	Magnesium (Mg ²⁺)		3750	mg/L		
	Calcium (Ca ²⁺)		7500	mg/L		