

## **Custom Solution - 11 components in HNO3 5%**

Art. ID LM24-CUS-72206

Unit 1 x 100 mL bottle

Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

### **Description**

Conservation 24 months. (12 months after opening)

| Text/Information | Analyte/Parameter | CAS number   | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|--------------|---------------------|------|--------|--------|
|                  | Antimony (Sb)     | [7440-36-0]  | 100                 | mg/L |        |        |
|                  | Arsenic (As)      | [7440-38-2]  | 100                 | mg/L |        |        |
|                  | Molybdenum (Mo)   | [7439-98-7]  | 100                 | mg/L |        |        |
|                  | Tin (Sn)          | [7440-31-5]  | 100                 | mg/L |        |        |
|                  | Tellurium (Te)    | [13494-80-9] | 100                 | mg/L |        |        |
|                  | Titanium (Ti)     | [7440-32-6]  | 100                 | mg/L |        |        |
|                  | Vanadium (V)      | [7440-62-2]  | 100                 | mg/L |        |        |
|                  | Thallium (Tl)     | [7440-28-0]  | 100                 | mg/L |        |        |
|                  | Silicon (Si)      | [7440-21-3]  | 1000                | mg/L |        |        |
|                  | Phosphorus (P)    | [7723-14-0]  | 1000                | mg/L |        |        |
|                  | Boron (B)         | [7440-42-8]  | 1000                | mg/L |        |        |