

Custom Solution - 24 components 5 mg/L each of in HNO3 5%

Art. ID LM24-CUS-72211

Unit 1 x 100 mL bottle

Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation 12 months. (6 months after opening)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	5	mg/L		
	Aluminium (Al)	[7429-90-5]	5	mg/L		
	Arsenic (As)	[7440-38-2]	5	mg/L		
	Boron (B)	[7440-42-8]	5	mg/L		
	Bismuth (Bi)	[7440-69-9]	5	mg/L		
	Calcium (Ca)	[7440-70-2]	5	mg/L		
	Cadmium (Cd)	[7440-43-9]	5	mg/L		
	Cobalt (Co)	[7440-48-4]	5	mg/L		
	Chromium (Cr)	[7440-47-3]	5	mg/L		
	Copper (Cu)	[7440-50-8]	5	mg/L		
	Iron (Fe)	[7439-89-6]	5	mg/L		
	Gallium (Ga)	[7440-55-3]	5	mg/L		
	Germanium (Ge)	[7440-56-4]	5	mg/L		
	Indium (In)	[7440-74-6]	5	mg/L		
	Magnesium (Mg)	[7439-95-4]	5	mg/L		
	Manganese (Mn)	[7439-96-5]	5	mg/L		
	Sodium (Na)	[7440-23-5]	5	mg/L		
	Nickel (Ni)	[7440-02-0]	5	mg/L		
	Lead (Pb)	[7439-92-1]	5	mg/L		
	Antimony (Sb)	[7440-36-0]	5	mg/L		
	Silicon (Si)	[7440-21-3]	5	mg/L		
	Tungsten (W)	[7440-33-7]	5	mg/L		
	Zinc (Zn)	[7440-66-6]	5	mg/L		
	Zirconium (Zr)	[7440-67-7]	5	mg/L		