

1-Benzhydryl-4-nitrosopiperazine Solution (1 mg/mL in Methanol)
(1-benzhydryl-4-nitrosopiperazine)

Art. ID USP-1A04130
Unit 1 mL
Deliverydetails De Minimis UN1230 (6.1) II E2

Description

Note: PAI products are different from official USP reference standards. These products are released using a separate process developed by USP's subject matter experts. The release process follows rigorous internal policies, standard operating procedures, and USP Quality Management System requirements. USP is an ISO9001:2015 certified facility.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	1-Benzhydryl-4-nitrosopiperazine	[1698-25-5]	1	mg/mL		