

### **Custom Solution - W 10 mg/L in HNO3 2%**

Art. ID LM24-CUS-72217  
Unit 1 x 100 mL bottle  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

#### Description

Conservation 12 months. (6 months after opening)

| Text/Information | Analyte/Parameter | CAS number  | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
|                  | Tungsten (W)      | [7440-33-7] | 10                  | mg/L |        |        |