

Custom Solution - 6 components 10 µg/L each of in HNO3 2%, HF 0.1%

Art. ID LM24-CUS-72218
Unit 1 x 100 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN2922 (6.1) III E1

Description

Conservation 8 months. (4 months after opening)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Germanium (Ge)	[7440-56-4]	10	µg/L		
	Zirconium (Zr)	[7440-67-7]	10	µg/L		
	Niobium (Nb)	[7440-03-1]	10	µg/L		
	Tungsten (W)	[7440-33-7]	10	µg/L		
	Tantalum (Ta)	[7440-25-7]	10	µg/L		
	Hafnium (Hf)	[7440-58-6]	10	µg/L		