

Custom Solution - Dicofol 100 µg/mL in Methanol

Art. ID LM24-CUS-71004-1X2ML
Unit 1 x 2 mL ampoule
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

Description

Conservation: 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Dicofol	[115-32-2]	100	µg/mL		