

Custom Solution - 10 components, in HNO3 2%

Art. ID LM24-CUS-71239-100ML
Unit 100 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Conservation: 12 months (originally closed bottle), (6 months after opening)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Phosphorus (P)	[7723-14-0]	1000	mg/L		
	Potassium (K)	[7440-09-7]	500	mg/L		
	Nickel (Ni)	[7440-02-0]	500	mg/L		
	Aluminium (Al)	[7429-90-5]	100	mg/L		
	Copper (Cu)	[7440-50-8]	100	mg/L		
	Manganese (Mn)	[7439-96-5]	100	mg/L		
	Barium (Ba)	[7440-39-3]	20	mg/L		
	Calcium (Ca)	[7440-70-2]	20	mg/L		
	Magnesium (Mg)	[7439-95-4]	20	mg/L		
	Zinc (Zn)	[7440-66-6]	20	mg/L		