

Soybean BASF Non-Modified Leaf DNA 0707-A9

Art. ID	AOCS-0707-A9
Unit	approx. 10 µg
Deliverydetails	No Dangerous Good /not restricted

Description

AOCS 0707-A9 was prepared from soybean genomic DNA extracted from clean non-modified soybean leaves provided by BASF. AOCS 0707-A9 is available in sterile, 0.5 ml skirted, screw-cap, self-sealing tubes containing approximately 10 µg extracted DNA. /// Certified Presence: FG72, A2704-12 and A5547-127 DNA absent /// Test Method: event-specific real-time PCR

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
FG72	FG72		<0,05	ng/µg	event-spec ific real- time PCR	
A2704-12	A2704-12		<0,05	ng/µg	event-spec ific real- time PCR	
A5547-127 DNA absent	A5547-127 DANN		<0,05	ng/µg	event-spec ific real- time PCR	