

Custom Solution - 24 components in HNO3 5%

Art. ID LM24-CUS-71799-5X100ML
Unit 5 x 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation 12 months (originally closed bottle). (6 months after opening)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Arsenic (As)	[7440-38-2]	50	mg/L		
	Silver (Ag)	[7440-22-4]	50	mg/L		
	Boron (B)	[7440-42-8]	50	mg/L		
	Barium (Ba)	[7440-39-3]	250	mg/L		
	Beryllium (Be)	[7440-41-7]	50	mg/L		
	Bismuth (Bi)	[7440-69-9]	50	mg/L		
	Cadmium (Cd)	[7440-43-9]	5	mg/L		
	Cerium (Ce)	[7440-45-1]	50	mg/L		
	Cobalt (Co)	[7440-48-4]	50	mg/L		
	Chromium (Cr)	[7440-47-3]	250	mg/L		
	Copper (Cu)	[7440-50-8]	250	mg/L		
	Gadolinium (Gd)	[7440-54-2]	2,5	mg/L		
	Lithium (Li)	[7439-93-2]	50	mg/L		
	Manganese (Mn)	[7439-96-5]	1000	mg/L		
	Molybdenum (Mo)	[7439-98-7]	50	mg/L		
	Nickel (Ni)	[7440-02-0]	250	mg/L		
	Lead (Pb)	[7439-92-1]	250	mg/L		
	Rubidium (Rb)	[7440-17-7]	50	mg/L		
	Selenium (Se)	[7782-49-2]	50	mg/L		
	Strontium (Sr)	[7440-24-6]	50	mg/L		
	Thallium (Tl)	[7440-28-0]	5	mg/L		
	Uranium (U)	[7440-61-1]	5	mg/L		
	Vanadium (V)	[7440-62-2]	50	mg/L		
	Zinc (Zn)	[7440-66-6]	250	mg/L		