

Custom Solution - 12 components in HNO3 2%

Art. ID LM24-CUS-70640-100ML
Unit 100 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Conservation 12 months. (6 months after opening)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	1	mg/L		
	Aluminium (Al)	[7429-90-5]	10	mg/L		
	Boron (B)	[7440-42-8]	10	mg/L		
	Calcium (Ca)	[7440-70-2]	1000	mg/L		
	Copper (Cu)	[7440-50-8]	10	mg/L		
	Iron (Fe)	[7439-89-6]	10	mg/L		
	Potassium (K)	[7440-09-7]	100	mg/L		
	Magnesium (Mg)	[7439-95-4]	200	mg/L		
	Sodium (Na)	[7440-23-5]	100	mg/L		
	Phosphorus (P)	[7723-14-0]	10	mg/L		
	Silicon (Si)	[7440-21-3]	50	mg/L		
	Zinc (Zn)	[7440-66-6]	50	mg/L		