

### **Custom Solution - Dibenzo(a,h)anthracene D14 100 µg/mL in Toluene**

Art. ID                      LM24-CUS-72239  
Unit                         3 x 1 mL ampoule  
Deliverydetails            De Minimis UN1294 II E2

#### Description

Conservation 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Dibenz(a,h)anthracene D 14	[13250-98-1]	100	µg/mL		