

Iso-octane

Art. ID LM24-OP-4008109-25ML
Unit 25 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1262 II E2

Description

Shelf-Life: 24 months /// Produced under ISO 17034 and ISO 17025

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|---------------------------------------|------------|---------------------|------|--------|--------|
| | 2,2,4-Trimethylpentane (Isooctane) | [540-84-1] | | | | |