

### **Custom Solution - Zn 2mg/L in HNO3 2%**

Art. ID LM24-CUS-72209-1X100ML  
Unit 1 x 100 mL bottle  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

#### **Description**

Conservation 12 months. (6 months after opening)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Zinc (Zn)	[7440-66-6]	2	mg/L		