

Titanium grade, CRM, Disc of Ø 5 mm

Art. ID	TL-3101
Unit	50 discs
Deliverydetails	No Dangerous Good /not restricted

Description

The sample TL-3101 is a titanium disc of 5 mm diameter and an average sample mass of 0,17g per pin. This material is packaged in a bottle of 50 discs.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Oxygen (O)	[7782-44-7]	0,1721 ± 0,0061	%		
	Nitrogen (N)	[7727-37-9]	0,0039 ± 0,0008	%		
	Hydrogen (H)	[1333-74-0]	0,00127 ± 0,0012	%		