

CHEMPLEX, XRF Sample Cups with TrimLess® Sleeve, Internal Overflow Reservoir with Vented Cap, Double Open Ends & Stackable, 1.24" (31.5mm) Dia X 1.15" (29.2mm) Height

Art. ID	CX-2132S
Unit	100 sets per package
Deliverydetails	No Dangerous Good /not restricted

Description

Snap-On Cap establishes a series of physical baffles that provide resistance to powdered sample withdrawal under evacuated conditions. Vent hole equalizes pressure differentials within the sample cup and sample chamber under helium flushes and vacuum operating conditions /// Features: TrimLess® sleeve completely envelops the thin-film/Double open-ended/Top sample loading/Internal Overflow Reservoir/Vented friction-fitting cap/He Purge support for P2