

A2704-12 soybean leaf tissue genomic DNA

Art. ID AOCS-0707-B15
Unit approx. 7.2 µg
Deliverydetails No Dangerous Good /not restricted

Description

AOCS-0707-B15 has been prepared by AOCS from homozygous genomic soybean DNA extracted from clean leaves. AOCS 0707-B15 is available in sterile, 0.5 mL skirted screw-cap self-sealing tubes containing approximately 7.2 µg of DNA. Each sample contains approximately 45.6 µL at a concentration of 150.7 ng/µL.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Homozygous A2704-12 DN A2704-12 DNA present			>999,99	ng/µg	Event specific real-time PCR	