

## **Azaspiracid-1 standard solution in Methanol**

Art. ID	NRC-CRM-AZA1-c
Unit	approx. 0,5 mL
Deliverydetails	De Minimis UN1230 (6.1) II E2

### Description

NRC-CRM-AZA1-c is a certified calibration solution for analytical method development and accurate quantitation of AZA1. The concentration is suitable for preparing a dilution series for calibration of instruments such as liquid chromatography with detection by mass spectrometry (LC-MS), as well as for spiking control samples for recovery experiments.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
15 - 30 °C	Azaspiracid-1	[214899-21-5]	1,44 ± 0,07	µmol/L		
15 - 30 °C	Azaspiracid-1	[214899-21-5]	1,22 ± 0,06	µg/mL		
	Azaspiracid-1	[214899-21-5]	1,54 ± 0,08	µg/g		