

Methyl Dithienyl glycolate (Methyl 2-hydroxy-2,2-di(thiophen-2-yl)acetate)

Art. ID USP-1A05100
Unit 25 mg

Description

Note: PAI products are different from official USP reference standards. These products are released using a separate process developed by USP's subject matter experts. The release process follows rigorous internal policies, standard operating procedures, and USP Quality Management System requirements. USP is an ISO9001:2015 certified facility.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Methyl Di(2-thienylglycolate)	[26447-85-8]				