

N-Nitroso Quinapril Solution (1 mg/mL in Acetonitrile)
((S)-2-(N-((S)-1-ethoxy-1-oxo-4-phenylbutan-2-yl)-N-nitroso-L-alanyl)-1,2,3,4-tetrahydroisoquinoline-3-carboxylic acid)

Art. ID USP-1A04180
Unit 1 mL
Deliverydetails De Minimis UN1648 II E2

Description

Note: PAI products are different from official USP reference standards. These products are released using a separate process developed by USP's subject matter experts. The release process follows rigorous internal policies, standard operating procedures, and USP Quality Management System requirements. USP is an ISO9001:2015 certified facility.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	N-Nitroso Quinapril	[n/a]	1	mg/mL		