

### **Custom Solution - 14 components 100 mg/L each of in HNO3 5%**

Art. ID LM24-CUS-72282  
Unit 250 mL bottle  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

#### Description

Conservation 24 months. (12 months after opening)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	100	mg/L		
	Arsenic (As)	[7440-38-2]	100	mg/L		
	Beryllium (Be)	[7440-41-7]	100	mg/L		
	Cadmium (Cd)	[7440-43-9]	100	mg/L		
	Cobalt (Co)	[7440-48-4]	100	mg/L		
	Chromium (Cr)	[7440-47-3]	100	mg/L		
	Copper (Cu)	[7440-50-8]	100	mg/L		
	Molybdenum (Mo)	[7439-98-7]	100	mg/L		
	Nickel (Ni)	[7440-02-0]	100	mg/L		
	Lead (Pb)	[7439-92-1]	100	mg/L		
	Antimony (Sb)	[7440-36-0]	100	mg/L		
	Selenium (Se)	[7782-49-2]	100	mg/L		
	Thallium (Tl)	[7440-28-0]	100	mg/L		
	Vanadium (V)	[7440-62-2]	100	mg/L		