

Dinophysistoxin-2 (DTX2) in Methanol

Art. ID	NRC-CRM-DTX2-c
Unit	0,5 mL
Deliverydetails	De Minimis UN1230 (6.1) II E2

Description

NRC-CRM DTX2-b is a certified calibration solution for the determination of dinophysistoxin-2 (DTX2) by chemical analysis methods. Each ampoule contains ~0.5 mL of a solution DTX2 in methanol at a concentration of 4.7 µmol/L.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Dinophysistoxin-2 (DTX2))	[139933-46-3]	4,86 ± 0,23	µg/g		
15 - 30 °C	Dinophysistoxin-2 (DTX2))	[139933-46-3]	3,83 ± 0,18	µg/mL		
15 - 30 °C	Dinophysistoxin-2 (DTX2))	[139933-46-3]	4,76 ± 0,23	µmol/L		