

Custom Solution - 20 components in HNO₃ 5%

Art. ID LM24-CUS-72292
Unit 1 x 100 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation 24 months(originally closed bottle), (12 months after opening).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Beryllium (Be)	[7440-41-7]	100	mg/L		
	Thallium (Tl)	[7440-28-0]	100	mg/L		
	Titanium (Ti)	[7440-32-6]	500	mg/L		
	Vanadium (V)	[7440-62-2]	500	mg/L		
	Chromium (Cr)	[7440-47-3]	500	mg/L		
	Manganese (Mn)	[7439-96-5]	500	mg/L		
	Cobalt (Co)	[7440-48-4]	500	mg/L		
	Nickel (Ni)	[7440-02-0]	500	mg/L		
	Strontium (Sr)	[7440-24-6]	500	mg/L		
	Molybdenum (Mo)	[7439-98-7]	500	mg/L		
	Tin (Sn)	[7440-31-5]	500	mg/L		
	Tellurium (Te)	[13494-80-9]	500	mg/L		
	Antimony (Sb)	[7440-36-0]	500	mg/L		
	Copper (Cu)	[7440-50-8]	500	mg/L		
	Arsenic (As)	[7440-38-2]	250	mg/L		
	Selenium (Se)	[7782-49-2]	250	mg/L		
	Lead (Pb)	[7439-92-1]	250	mg/L		
	Silver (Ag)	[7440-22-4]	250	mg/L		
	Cadmium (Cd)	[7440-43-9]	100	mg/L		
	Uranium (U)	[7440-61-1]	100	mg/L		