

Custom Solution - 10 components in HNO3 5%

Art. ID	LM24-CUS-70897-100ML
Unit	1 x 100 mL bottle
Deliverydetails	ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation: 24 months (originally closed bottle), 12 months after opening.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	1000	mg/L		
	Iron (Fe)	[7439-89-6]	1000	mg/L		
	Zinc (Zn)	[7440-66-6]	1000	mg/L		
	Barium (Ba)	[7440-39-3]	1000	mg/L		
	Boron (B)	[7440-42-8]	2500	mg/L		
	Sodium (Na)	[7440-23-5]	10000	mg/L		
	Phosphorus (P)	[7723-14-0]	5000	mg/L		
	Magnesium (Mg)	[7439-95-4]	10000	mg/L		
	Calcium (Ca)	[7440-70-2]	10000	mg/L		
	Potassium (K)	[7440-09-7]	10000	mg/L		