

## **Custom Solution - 7 components in Water**

Art. ID	LM24-CUS-71338-50ML
Unit	50 mL bottle
Deliverydetails	No Dangerous Good /not restricted

### Description

Conservation 24 months (originally closed bottle) - (12 months after opening).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Nitrate (NO <sub>3</sub> -)		2500	mg/L		
	Nitrite (NO <sub>2</sub> -)		250	mg/L		
	Chloride (Cl-)		10000	mg/L		
	Sulfate (SO <sub>4</sub> <sup>2-</sup> )		5000	mg/L		
	Fluoride (F-)		500	mg/L		
	Bromide (Br-)		250	mg/L		
	Phosphate (PO <sub>4</sub> <sup>3-</sup> )		500	mg/L		