

Custom Solution - 4 components in HNO3 5%

Art. ID LM24-CUS-72303
Unit 1 x 100 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation 12 months (originally closed bottle), (6 months after opening).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	30	mg/L		
	Manganese (Mn)	[7439-96-5]	30	mg/L		
	Iron (Fe)	[7439-89-6]	30	mg/L		
	Zinc (Zn)	[7440-66-6]	10	mg/L		