

Custom Solution - 16 components in HNO3 5%

Art. ID	LM24-CUS-72315
Unit	1 x 500 mL bottle
Deliverydetails	ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation 12 months (originally closed bottle), (6 months after opening).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	50	mg/L		
	Aluminium (Al)	[7429-90-5]	50	mg/L		
	Boron (B)	[7440-42-8]	50	mg/L		
	Barium (Ba)	[7440-39-3]	50	mg/L		
	Calcium (Ca)	[7440-70-2]	1000	mg/L		
	Cadmium (Cd)	[7440-43-9]	50	mg/L		
	Chromium (Cr)	[7440-47-3]	50	mg/L		
	Copper (Cu)	[7440-50-8]	100	mg/L		
	Iron (Fe)	[7439-89-6]	100	mg/L		
	Potassium (K)	[7440-09-7]	100	mg/L		
	Magnesium (Mg)	[7439-95-4]	200	mg/L		
	Manganese (Mn)	[7439-96-5]	50	mg/L		
	Sodium (Na)	[7440-23-5]	500	mg/L		
	Nickel (Ni)	[7440-02-0]	50	mg/L		
	Lead (Pb)	[7439-92-1]	50	mg/L		
	Zinc (Zn)	[7440-66-6]	100	mg/L		