

Custom Solution - 8 components in HNO3 5%

Art. ID LM24-CUS-72319
Unit 1 x 100 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation 24 months. (12 months after opening)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	500	mg/L		
	Barium (Ba)	[7440-39-3]	500	mg/L		
	Boron (B)	[7440-42-8]	1000	mg/L		
	Iron (Fe)	[7439-89-6]	1000	mg/L		
	Phosphorus (P)	[7723-14-0]	1000	mg/L		
	Strontium (Sr)	[7440-24-6]	1000	mg/L		
	Zinc (Zn)	[7440-66-6]	1000	mg/L		
	Potassium (K)	[7440-09-7]	10000	mg/L		