

Custom Solution - 18 components 100 mg/L each of in Acetonitrile

Art. ID LM24-CUS-72320
Unit 3 x 1 mL ampoule
Deliverydetails De Minimis UN1648 II E2

Description

Conservation 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2-Phenylphenol	[90-43-7]	100	mg/L		
	Benzovindiflupyr	[1072957-71-1]	100	mg/L		
	Binapacryl	[485-31-4]	100	mg/L		
	Bixafen	[581809-46-3]	100	mg/L		
	Diniconazole	[83657-24-3]	100	mg/L		
	Diphenylamine	[122-39-4]	100	mg/L		
	Fenvalerate	[51630-58-1]	100	mg/L		
	Fipronil-sulfone	[120068-36-2]	100	mg/L		
	Fluazifop-methyl	[69335-90-6]	100	mg/L		
	Fluoxastrobin	[361377-29-9]	100	mg/L		
	Flurprimidol	[56425-91-3]	100	mg/L		
	Isofetamid	[875915-78-9]	100	mg/L		
	Pyraflufen (free acid)	[129630-17-7]	100	mg/L		
	S-Metolachlor	[87392-12-9]	100	mg/L		
	Tiocarbazil	[36756-79-3]	100	mg/L		
	Dichlofluanid	[1085-98-9]	100	mg/L		
	Parathion-ethyl	[56-38-2]	100	mg/L		
	Fenitrothion	[122-14-5]	100	mg/L		