

Custom Solution - 2 components 100 µg/mL each of in Methanol

Art. ID LM24-CUS-72321
Unit 3 x 1 mL ampoule
Deliverydetails De Minimis UN1230 (6.1) II E2

Description

Conservation 12 months (in originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Triphenylmethan	[519-73-3]	100	µg/mL		
	Triphenylphosphate	[115-86-6]	100	µg/mL		